

COVID-19 Screening (Express Service) (CS1) (by Virtual Care instruction)

Last update on: 27-7-2020

Reserve time slot and purchase



Attend Virtual Care instruction



Collect the deep throat saliva collection kit at
designated UMP medical centre



Deep throat saliva collection tutorial



Collect deep throat saliva specimen

Time of specimen collection:
After getting up and before brushing teeth,
eating or drinking



- Collect and submit specimen on the same day*
- Specimen submission arrangement: specimen can be submitted to designated UMP medical centres*
- Under general circumstances, the result will be emailed within 24 hours at the earliest.*

Enquiry Hotline : 3110 2785

- Services provided on a first-come-first-served basis.
- UMP reserves the right of final decision and interpretation.
- *Please refer to the service booking page for details. (<https://vcbooking.ump.com.hk/v2/#book>)